Little Doctor LD-152U

EN ULTRASONIC NEBULIZER LD. Instruction Manual.

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
















1. If you are not using the device, disconnect it from the mains by pulling out the plug from the socket.

- Observe the safety regulations for electrical appliances, in particular:
 - Never immerse the device in water;
 - Install the device vertically on a flat, stable surface;
 - Position the device so as not to block the ventilation holes on the bottom of the device;
 - Do not use this device in explosive environments, e.g., near flammable anesthetics or inside an oxygen chamber;
 - Do not touch the power plug with wet hands;
 - Do not place the device in water, under a water drain or in a stall shower. Do not use when bathing;
 - Do not touch the device if it has fallen into water. Disconnect it from the mains immediately;
 - The device and especially its accessories should be kept out of the reach of children, as they contain small parts that can be swallowed;
 - Time of contact with the power cord to remove the plug from the socket;
 - The power cord must not be in contact with hot or heated surfaces;
 - Store and use the device indoors, away from heat sources. After use, it is recommended to store the device away from dust and direct sunlight.

- If the device does not work, see the section "TROUBLESHOOTING ALGORITHM AND REMEDIES".
- Do not insert foreign objects into the openings of the device.
- The device is not intended for outdoor use.
- None of the electrical or mechanical parts are intended to be repaired by end users. Do not open the device yourself! For repairs, contact exclusively the authorized service centers listed in this Manual. Therefore, it should be used as an aerosol therapy system. Any use other than that for which the device is intended is considered inappropriate and manufacturer cannot be held responsible for damage caused by incorrect, erroneous and/or unreasonable use, or if the device is used in electrical systems that do not comply with current regulations;
- None of the electrical or mechanical parts are intended to be repaired by end users. Do not open the device yourself! For repairs, contact exclusively the authorized service centers listed in this Manual. Therefore, it should be used as an aerosol therapy system. Any use other than that for which the device is intended is considered inappropriate and manufacturer cannot be held responsible for damage caused by incorrect, erroneous and/or unreasonable use, or if the device is used in electrical systems that do not comply with current regulations;
- The instrument requires special precautions with regard to electromagnetic compatibility and must be installed and used in accordance with the information in the "ELECTROMAGNETIC COMPATIBILITY INFORMATION" section.
- Keep the device accessories out of the reach of children. Children and persons with learning difficulties or lack of the necessary skills should only use the device under the strict supervision of fully capable adults who have carefully read this Manual.
- Use only the accessories designed for the LD-152U and described in this Manual.
- Parts in contact with patients such as masks, mouthpieces or nosepieces comply with ISO 10993-1 and do not cause allergic reactions or skin irritation.
- The device and its parts are biocompatible in accordance with the requirements of the EN 60601-1 standard.
- Do not use the device, no additional instructions are required other than those indicated in this Manual.
- Do not use the device in conditions other than those specified in this Manual may seriously compromise the safety of the device and its technical characteristics;
- The parts of the device that are in contact with drugs are made of high stable chemical-resistant thermoplastic polymers. Such materials have been tested with drugs commonly used for nebulizer therapy (salbutamol, beclomethasone dipropionate, acetylcysteine, budesonide, ambroxol) and no effect was noted on the materials. However, such an effect cannot be completely ruled out, given the variety and constant evolution of the drugs used. Remember, such an effect can occur even if the drug is used as soon as possible after opening the package;
- Avoid long-term storage of the drug in the container and clean it immediately after each use;
- Do not use drugs in damaged packaging;
- Do not forget:
 - The nature of inhalation (through the mouth using a mouthpiece or through the nose using a mask), timing (maximum 10-15 minutes as usual), and frequency, as well as the nebulizer solutions used, should be prescribed by YOUR DOCTOR.
 - Inhalation with a mask is used in children who cannot use the mouthpiece.
 - Solutions for inhalation should be prepared under sterile conditions using 0.9% sodium chloride as a solvent. You should not use tap water, even after boiling. The vessels in which the solution is prepared are previously disinfected by boiling.

⚠ Attention! This device is not intended for inhalation anesthesia or mechanical ventilation.

OPERATING PROCEDURES

⚠ Important! Before using the device for the first time, it must be thoroughly cleaned, as described in paragraph 1 of the "Maintenance, storage, repair, and disposal" section

- Remove the nebulizer chamber from the main unit of the device, holding it by the tabs on the nebulizer chamber cap (Figure 2a)
- Fill the container with water exactly up to the "WATER LEVEL" mark on the level sensor. If the water level is less than the permissible level, the sensor will not allow the piezoelectric element to be turned on, while the red indicator will be on «**△**».
- Place the nebulizer chamber vertically on a table. Detail the top of the nebulizer chamber from the chamber. To pull the cap of the nebulizer chamber off counterclockwise, so that the marker «**▶▶▶**» points to «**⬇**» (Figure 2b). Pull the cap of the nebulizer chamber up and remove it.
- Fill the container with no more than 12 mL and at least 2 mL of a nebulizer solution.

⚠ Attention! Do not turn on the device with an empty drug container. Do not add more than 12 mL of nebulizer solution into the container.

- Close the nebulizer chamber cap back by turning it clockwise so that the marker «**▶▶▶**» points to «**⬆**» (Figure 2c).
- Place the nebulizer chamber back into the main unit of the device.
- Connect the tube with a mask or mouthpiece to the nebulizer chamber. If necessary, you can lengthen the tube by stretching it or connecting the two tubes with a sleeve.

It should be borne in mind that increasing the length of the inhalation tubes leads to an increase in the loss of the nebulizer solution.

- ⚠ Important!** Each patient is advised to use an individual mouthpiece and/or mask.
- Set the controls to their original positions:
 - regulator «**⊗**» to position «0»;
 - airflow regulator «**⊕**» to position «1»;
 - air flow power plug into a socket.
- ⚠ Important!** Nearby the device. Put on a mask or put a mouthpiece in your mouth. DEVICE IS READY FOR INHALATION.

⚠ Attention! For effective inhalation, you need to calm down, relax, and sit up straight.

If you breathe in bed, place a pillow under your back so that your back is as straight as possible. An awkward posture and uneven inhalation can block some parts of the lungs and the airway.

- To start inhalation, turn the regulator «**⊗**» clockwise to the position "1", the green indicator «**⊕**» will light up, and the device will begin to produce an aerosol.
- In the initial period of inhalation, adjust the optimal level of aerosol formation intensity using the regulator «**⊗**». The intensity of aerosol formation, and consequently, the nebulizer solution flow rate depends on the viscosity of the solution used. For example, when using 0.9% sodium chloride as a nebulizer solution, the dependence of the nebulizer solution flow time on the regulator position «**⊗**» is shown in the table.

⚠ Attention! Do not turn on the device with an empty drug container. Do not add more than 12 mL of nebulizer solution into the container.

⚠ Important! Each patient is advised to use an individual mouthpiece and/or mask.

NOTE 1: If you are using a nebulizer for children, too much aerosol flow can make it difficult for the baby to breathe. In this case, it is recommended to reduce the aerosol flow with the air flow regulator, which will increase the inhalation time.

NOTE 2: During inhalation, condensation may form under the cap of the device.

- After using:
 - a) turn off the device by setting the regulator «**⊗**» to the position "0";
 - b) remove the power plug from the socket;
 - c) disconnect the nebulizer chamber from the main unit, remove the upper cap of the nebulizer chamber and drain the remaining nebulizer solution;
 - d) clean the water from the container;
 - e) wash the device with a soft cloth and let it dry;
 - f) move the device to a dry place out of reach of children.

REPLACING THE CONTAINER OF SOLUTION

Inspect the container of the solution before each use. In case of violation of integrity, deformation or other damage, replace the container with a new one.

Proceed as follows to replace the container:

1. Remove the nebulizer chamber from the main unit (pos. 1 and 2 of the "PARTS AND COMPONENTS" section).
2. Unscrew the lower part of the nebulizer chamber (pos. 1.6).
3. First remove the gasket (pos. 1.4) and then the container (pos. 1.5) from the bottom of the nebulizer chamber.
4. Place the new container in the bottom of the nebulizer chamber and install the gasket on top. Make sure the container and gasket fit tight against the rails.
5. Screw the lower part of the nebulizer chamber to the middle part of the chamber (pos. 1.3).
6. Install the nebulizer chamber into the main unit.

MAINTENANCE, STORAGE, REPAIR, AND DISPOSAL

1. Clean the device and all accessories regularly. Disinfect device accessories after each use. All device accessories are recommended to be treated by wiping, sprinkling, immersing, and soaking:

- a) In medical institutions, using special disinfectants. Mouthpieces can be boiled for 10 minutes or autoclaved up to 150 °C. After treatment, wipe all parts of the device dry with a soft cloth.
- b) At home, using a 3-6% hydrogen peroxide solution. The disinfection time is 30 min. The initial temperature of the solution is 50 °C. The time of the working solution is 1 day. Or using a 3-6% solution of table vinegar. The solution is prepared from 1 part of 3-6% table vinegar and 3 parts of distilled water. Make enough solution to immerse the parts. The disinfection time is 30 min. Do not use working solution! After treatment, rinse all the device accessories under an abundant stream of hot water. Air-dry the accessories or blot with a soft, clean, lint-free cloth.

2. Do not expose air tubes and masks to boiling, autoclaving, or dishwasher cleaning!
3. The device must be protected from direct sunlight and shock.
4. Do not store or use the device in the vicinity of heaters and open flames.
5. Protect the device from contamination.
6. Do not allow the device to come into contact with aggressive solutions.
7. If necessary, carry out repairs only with specialized contractors.
8. The life cycle of the drug container is 6 months from the beginning of the operation. The manufacturer has not established special conditions for disposal.

WARRANTY OBLIGATIONS

This device has a warranty period of 36 months from the date of sale. The warranty does not cover consumables (nebulizer solution containers, masks, mouthpieces, tubes, etc.).

The manufacturer is not responsible for any accidental or collateral damage resulting from actions to modify or repair the device, performed not in authorized service centers, or if any of the components of the device has been damaged as a result of an accident or misuse. Any unprofessional modification/repair of the device will void the warranty, and the manufacturer does not guarantee further compliance of the device with the requirements of the Medical Device Directive 93/42/EEC (and subsequent changes) and its regulatory documents.

Warranty obligations shall be formalized with a warranty card when the device is sold to a customer. The addresses of the organizations providing warranty service are indicated on the website www.littledoctor.pl.

⚠ Attention! It is not allowed to use nebulizer solutions containing ether, oils, or suspended particles (suspensions), including decoctions and herbal infusions. All types of standard nebulizer solutions in liquid form for nebulizer therapy produced by pharmaceutical companies are recommended for use.

⚠ Attention! This device is not intended for inhalation with solutions of hormonal drugs and antibiotics, because the work of the piezoelectric element leads to their destruction.

1. After opening the package, check the integrity of the device, paying particular attention to the presence of damage to plastic parts, which can make accessible parts of the under-voltage device, as well as damage to the power cable. Switch on the device only after making sure that there is no such damage. Perform this check before each use.
2. Before connecting the device, make sure that the electrical characteristics of the device correspond to those of the mains.

TECHNICAL SPECIFICATIONS

LD-152U	ULTRASONIC
Power consumption, no more	40 VA
Aerosol production performance	2 ml/min *
Mean diameter of aerosol particles (MMAD)	3 μm*
Ultrasonic frequency	1.7 MHz ± 10%
Water temperature	< 60 °C
Volume of the nebulizer solution container	12 mL
Residual volume of nebulizer solution, no more	1 mL
Time of continuous use	50 min
Cooldown time of the device	at least 40 min
Power supply:	~230V 50Hz
Protection against electric shock	Product type BF
Device operational conditions:	
Ambient temperature	10 °C to 40 °C
Humidity	15% to 85% Rh
Storage and transportation conditions for the device:	
Ambient temperature	-30 °C to 40 °C
Humidity	15% to 85% Rh
Weight of the device (without packaging), no more	1,515 kg
Overall dimensions of the electronic unit	250x104x160 mm
Year and month of production	Indicated on the body of the device in the serial number as «YYMMSSZZXXX», where YY is the year, and MM is the month of manufacture.

* - data obtained by Little Doctor International (S) Pte. Ltd.

SYMBOL EXPLANATION

- ⚠ Important: Read the Operation Manual
- CE 0123 Compliance with Directive 93/42/EEC
- Complies with Customs Union technical regulations
- Product not sterile
- Serial number
- Protection class II
- Product type BF
- Manufacturer
- Conditions for storage, transportation, and use
- Representative in the European Union
- Aerosol formation intensity regulator
- Aerosol formation intensity regulator
- Airflow regulator
- Indicator of insufficient water level in the container
- ON/OFF indicator of power supply

The revision date of this Operation Manual is shown on the last page as E00X/YYMM/NN, where YY is the year, MM is the revision month, and NN is the revision number.

TROUBLESHOOTING ALGORITHM AND REMEDIES

In the absence of aerosol formation, troubleshooting is recommended according to the following order:

1	Indicator is off	Indicator is on
1	Check that the electrical plug is firmly connected to the socket. 2. Check the position of the regulator « ⊗ ». The regulator must not be in the «OFF» position 3. Check fuses and replace if necessary If the above steps do not correct the problem, contact a SPECIALIZED REPAIR CONTRACTOR	1. Check the float of the water level sensor 2. Fill the container with water to the «WATER LEVEL» mark 3. Go to step 3
2	Check the volume of the nebulizer solution in the nebulizer chamber (more than 1 mL).	There is sufficient nebulizer solution Go to step 4
3	Check the composition of the drug solution.	The nebulizer solution contains one of the listed components: ether, oils, including decoctions and herbal infusions. The nebulizer solution does NOT contain ether, oils, or suspended particles (suspensions) including decoctions and herbal infusions. It is not allowed to use nebulizer solutions containing ether, oils, or suspended particles (suspensions) including decoctions and herbal infusions. All types of standard nebulizer solutions in liquid form for nebulizer therapy produced by pharmaceutical companies are recommended for use. Contact a SPECIALIZED REPAIR CONTRACTOR.

⚠ Attention! At the same time, no more than 1 nebulizer solution should be installed in the nebulizer chamber

This device manufacturing is certified according to international standard ISO 13485. Device comply with the requirements of European Directive MDD 93/42/EEC.

Manufactured under control and for Little Doctor International (S) Pte. Ltd. (7500A BEACH ROAD #11-313 THE PLAZA SINGAPORE 199591. Postal address: Yishun Central PD, Box 9293 Singapore 917699).

Manufacturer: Little Doctor Electronic (Nantong) Co. Ltd. (No.8, Tongxing Road Economic & Technical Development Area, 226010 Nantong, Jiangsu, People's Republic of China).

Distributor in Europe: Little Doctor Europe Sp. z o.o., 576 Zawila Street, 30-390, Krakow, Poland, Sales Office phone: +48 12 2684746, 12 2684747, fax: +48 12 268 47 53. E-mail: biuro@littledoctor.pl, www.LittleDoctor.com

Authorized Representative in the EU: Little Doctor Europe Sp. z o.o., 576 Zawila Street Krakow 30-390 Poland.

For more information please contact www.littledoctor.pl

* - Measured by Little Doctor International (S) Pte. Ltd.

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